


[Loctite® Cured-In-Place Silicone Gasketing Properties Chart]

LOCTITE® PRODUCT	ITEM NUMBER	PACKAGE TYPE/SIZE	COLOR	EXTRUSION RATE (g/min.)	CURE METHOD	CURE SPEED*	HARDNESS	% ELONGATION	TENSILE STRENGTH (psi)	AGENCY APPROVALS
5039™ NUVA-SIL®	40438 40439	300 ml cartridge 40 lb. pail	Translucent/ Colorless/Fluorescent	222	UV or Visible/ Moisture	30 sec. @ 120 mW/cm ²	37	210	290	N/A
5050™ FASTGASKET® NEW	1212166 1212165	300 ml cartridge 40 lb. pail	Translucent	500	Ultraviolet	Ultraviolet	41	260	510	N/A
5950™ FASTGASKET® 	29287 18495	300 ml cartridge 40 lb. pail	Black	350	UV/Moisture	30 sec. @ 120 mW/cm ²	30	200	150	N/A
5951™ FASTGASKET®	18198	40 lb. pail	Clear	350	UV/Moisture	30 sec. @ 120 mW/cm ²	30	200	150	N/A
5960™ FASTGASKET®	929334	40 lb. pail	Clear Clear Blue	350	UV	30 sec. @ 120 mW/cm ²	50	260	600	N/A
5963™ PROCURE®	34337	50 lb. pail	Grey	250	Heat	10 min. @ 150°C	50	360	620	N/A
5964™ PROCURE®	34348 34347	300 ml cartridge 50 lb. pail	Brown	120	Heat/ Microwave	10 min. @ 150°C VFM 2 to 4 min.	32	530	635	N/A
3407™ OPTILOC®	35173	30 ml syringe	Black	5	Heat	60 min. @ 150°C	33	>600	>493	N/A

CURED-IN-PLACE

 = Indicates Worldwide Availability